

AVIDD ASAP: A71EV2A trunc08 A71EV2A-c012 Expression and Purification

PAGE23-00246

Author: Fairhead, Michael
Date Started: 2023-Feb-02
Experiment Started:
Projects: Expression;Purification;ASAP
Related Pages:
Referenced by:
Tags:

Title missing - double click to edit	
A71EV2A-c012/A71EV2A-k003/A71EV2A-e003/A71EV2A-p003	
A71EV2A _F_trunc08	A71EV2A-c012
MHHHHHHGSGDQEAKPSTEDLGDKKEGEYIKLKVIGQDSSEIHFKVKMTTHLKKLKESYCQRQGVPMNSLRFLFEGQRIADNHTPKELGMEEDVIEVYQEQTGGSGAIYVGNRVRVNRHLATHNDWANLVWEDSSRDLLVSSTTAQGCDTIARDCDQTGVYYCSSRRKHYPVSFSKPSLIFVEASEYYPARYQSHLMLAVGHSEPGDCGGILRCQHGVVGIVSTGGNGLVGFADVRDLLWLDEEAMEQ-	
Number of amino acids: 249	
Molecular weight: 27927.21	
Theoretical pI: 5.81	
Cut SUMO tag	
MHHHHHHGSGDQEAKPSTEDLGDKKEGEYIKLKVIGQDSSEIHFKVKMTTHLKKLKESYCQRQGVPMNSLRFLFEGQRIADNHTPKELGMEEDVIEVYQEQTGGSGAIYVGNRVRVNRHLATHNDWANLVWEDSSRDLLVSSTTAQGCDTIARDCDQTGVYYCSSRRKHYPVSFSKPSLIFVEASEYYPARYQSHLMLAVGHSEPGDCGGILRCQHGVVGIVSTGGNGLVGFADVRDLLWLDEEAMEQ-	
Number of amino acids: 105	
Molecular weight: 12069.48	
Theoretical pI: 5.96	
A71EV2A-c012	
SGAIYVGNRVRVNRHLATHNDWANLVWEDSSRDLLVSSTTAQGCDTIARDCDQTGVYYCSSRRKHYPVSFSKPSLIFVEASEYYPARYQSHLMLAVGHSEPGDCGGILRCQHGVVGIVSTGGNGLVGFADVRDLLWLDEEAMEQ-	
Number of amino acids: 144	
Molecular weight: 15875.74	
Theoretical pI: 5.58	
Extinction coefficients:	
Extinction coefficients are in units of M ⁻¹ cm ⁻¹ , at 280 nm measured in water.	
Ext. coefficient 28795	
Abs 0.1% (=1 g/l) 1.814, assuming all pairs of Cys residues form cystines	

Title missing - double click to edit

Expression

Transformed BL21[DE3]RR with **A71EV2A-c012** (/kanR)
Grew 100 mL o/n in SOC + Kan
Used 10 mL to inoculate 1 L FORM-TB + Kan (6 L total)
Grew 4h 37C 180 rpm shaking (OD600 = 1)
Grew 18C 180 rpm shaking 1h
Added 0.5 mM IPTG final conc.
Grew o/n 18C 180 rpm shaking
Harvested 4000g 12C 20 minutes
Froze pellet -80C (final wcv = 100 g/L 17 g total)

Formedium TB custom

12 g/L tryptone
24 g/L yeast extract
3.3 g/L Ammonium Sulphate (NH₄)₂SO₄
6.8g /L Potassium Dihydrogen Phosphate KH₂PO₄
7.1 g/L DiSodium Hydrogen Phosphate Na₂HPO₄ 7.1
0.15 g/L Magnesium Sulphate MgSO₄
0.03 g/L Trace Elements
55.85 g/L in MilliQ water added 20 mL 50 % glycerol
Autoclave
Add 1 mL 10% Antifoam 204 and Antibiotic

Title missing - double click to edit

Purification

Dissolve pellet in 4 mL per g using Lysis buffer (10 mM HEPES pH 7.5, 500 mM NaCl, 5 % glycerol, 0.5 mM TCEP, 1 % TX-100, 30 mM Imidazole, 0.5 mg/mL lysozyme, 0.01 mg/mL benzonase)
After dissolveing (1h RT stirring) incubate on ice 1h and then centrifuge 30000 g 1h 4C
Pass over 20 mL Ni-Sepharose-FF
Wash 3 x 100 mL Wash buffer (10 mM HEPES pH 7.5, 500 mM NaCl, 5 % glycerol, 0.5 mM TCEP, 30 mM Imidazole)
Elute 3 x 25 mL Elution Buffer (10 mM HEPES pH 7.5, 500 mM NaCl, 5 % glycerol, 0.5 mM TCEP, 500 mM Imidazole)

Pool peak fractions

50 mL at A280 of 12.7

Add 1:100 sumo protease (Senp1) and dialyse o/n against 1 L Wash Buffer

Pass over 20 mL Ni-Sepharose-FF (collect flow through)

Wash with 2 x 30 mL Wash Buffer 30 mM (collect flow through)

Wash with 2 x 30 mL Wash Buffer 60 mM (collect flow through)

Pool peak fractions

protein comes off predominately at 60 mM Imidazole but a lot of SUMO still present

Pass over 20 mL Ni-Sepharose-FF (collect flow through)

Wash with 2 x 30 mL Wash Buffer 60 mM (collect flow through)

Calculate yield and check on gel still some SUMO (similar MW and pI)

Concentrate to 9 mg/mL using 10,000 MWCO concentrator, struggles to get above A280 of 20 so lets not push it

Yield at this stage is approx 8 mL with A280 of 16, about 72 mg

Do SEC using 5 mL aliquots (approx. 40 mg) on 125 mL superose 12 pg using 10 mM HEPES pH 7.5, 500 mM NaCl, 5 % glycerol, 0.5 mM TCEP as mobile phase

Concentrate peak fractions to 0.79 mM (12.5 mg/mL, A280 of 22.5)

Do MS to confirm MW and flash freeze in LN2 and store in 100 µL single use aliquots (2.4 mL 30 mg total)

MS looks as expected very small signal for SUMO2 tag

A71EV2A-c012


Expected = 15875.74


Observed = 15876.12

SUMO-tag


Expected = 12069.48

Observed = 12069 but 6.8×10^2 compared to 0.9×10^7 for EVA712A





SDS-PAGE SEC FRACTIONS



Z:\Agilent_SGC_QTOF\Rod\230210_MS424H

